

Distal Torus Fracture in Forearm Pathway Splint Guide

How to size and apply a forearm splint:

1. Open all hook and loop closures, exposing the inside of the splint
2. Place the hand in the splint and slide the thumb through the thumb hole
3. Position splint with palmar side of the splint below the palmar crease.
 - a. Determine if stays (hard parts on top and bottom of splint) offer proper support and positioning. Remove splint if wrist position needs to be adjusted
4. Close the hook and loop straps to get a secure fit
 - a. Ensure the stays are directly over each other and offer support
 - b. Ensure the splint is tight enough to offer support without compromising circulation
 - c. If patient can freely move wrist in splint, it is not tight enough or is too big
5. Comfortably secure the wrist straps through the closures to hold splint in proper position
6. Trim excess strapping where needed.

Ordering:

Emergency Department	Richard/Velcro splint Capture in Oracle 6 inch: Right: 114950 Left: 114951 8 inch: R: 114980 L: 114981 10 inch: R: 114978 L: 114979
Urgent Care	Universal Wrist Splint Capture in Oracle #114981 #114979 #114950

Updated: 7/2020

DISTAL TORUS FRACTURE

SPLINT GUIDE



Children's Physicians order from CHHC

Splint Wrist Xxs Left Rj Ind.	5558-PWR-01-XXSL
Splint Wrist Xsm Left Rj Ind.	5558-PWR-01-XSL
Splint Wrist Left Pedi 8 Long	QTX507312
Splint Wrist Left X-Sm 8 Long	QTX507306
Splint Wrist Left Sm 8 Long	QTX507307
Splint Wrist Left Med 8 Long	QTX507308
Splint Wrist Left Lg 8 Long	QTX507309
Splint Wrist Left X-Lg 8 Long	QTX507310
Splint Wrist Xxs Right Rj Ind.	5558-PWR-01-XXSR
Splint Wrist Xsm Rt Rj Ind.	5558-PWR-01-XSR
Splint Wrist Rt Pedi 8 Long	QTX507311
Splint Wrist Rt X-Sm 8 Long	QTX507301
Splint Wrist Rt Sm 8 Long	QTX507302
Splint Wrist Rt Med 8 Long	QTX507303
Splint Wrist Rt Lg 8 Long	QTX507304
Splint Wrist Rt X-Lg 8 Long	QTX507305

Updated: 7/2020